



**Display at eye-level**

The display is positioned at eye-level to make it easier to read and simplify operation.



**Optimized air guide plate**

The optimized air guide plate reduces temperature deviations inside the system to a minimum.

# LABO 340/468

## PRO-ACTIVE

330 l /  
460 l

+2 to  
+20 °C

USB



### PRO-ACTIVE- Control:

Permanent, proactive monitoring of the performance data and alerts in the case of deviations – so that you can take countermeasures in good time before a fault threatens your chilled goods; world's most accurate temperature control in refrigerated areas thanks to two standard PT-1000 sensors.

■ **External housing** made from galvanised sheet steel (rust-proof), with white, anti-scratch powder coating. Length of the plug cable: approx. 2.8 m.

■ **Adjustable feet** to compensate uneven floors.

■ **Interior** made from smooth aluminium with colourless protective coating. Gibs (every 15 mm) for flexible interior arrangement.

■ **Interior** of LABO 340 PRO-ACTIVE consists of one black drawer and four plastic-coated shelves on overlays.

Interior of LABO 468 PRO-ACTIVE consists of four plastic-coated shelves on overlays.

■ **Extra-thick 55 mm energy saving insulation**, made from high-quality, compression-moulded and environmentally friendly material.

■ **Self-closing door** with easy-to-replace plastic magnetic seal frame, lockable.

■ **Door hinge** on the right-hand side (by default, see illustration), or the left-hand side, can be retrofitted.

■ **Forced-air cooling**, switches off automatically when you open the door, ensures a uniform temperature and minimizes temperature deviation.

■ **Automatic defrosting** with time limit and temperature monitoring.

■ **Condensate evaporates** in the refrigerating machine interior.

■ **Key switch** protected control panel against tampering. Membrane keyboard with digital temperature display and minimum/maximum temperature memory.

■ **Warning functions** with visual and audible alarm signal in the case of temperature deviations or other malfunctions. Door open alarm after 60 seconds.

■ **Alarms** can be forwarded using a potential-free contact (e.g., to a mobile phone with optional GSM-Module or to a control centre).

■ **Power failure alarm** (visual and audible alarm), the monitoring unit remains in operation for approx. 30 hours on battery power.

■ **Data documentation** can be read via USB interface with KIRSCH Datenet software.

■ **Antifreeze** against sub-zero temperatures.

■ **Statically ventilated** (LABO 340 PRO-ACTIVE) or **ventilation-enforced** (LABO 468 PRO-ACTIVE) refrigerating machine, hermetically sealed, energy saving, low-noise, easy to service.

## Specification

	LABO 340 PRO-ACTIVE	LABO 468 PRO-ACTIVE
Capacity	330 litres	460 litres
Temperature setting	approx. +0 °C to +15 °C	approx. +0 °C to +15 °C
Voltage	220 – 240 V, 50/60 Hz	220 – 240 V, 50/60 Hz
Power consumption	110 watts	170 watts
Normal consumption	0.65 kWh/24 h	1.00 kWh/24 h
Admissible ambient temperature	from +10 °C to +38 °C	from +10 °C to +38 °C
Heat emission (max.)	166 watts	280 watts
Exterior dimensions (including door handle and distance from wall)	w x d x h = 67 x 66 x 189 cm	w x d x h = 74 x 77,5 x 189 cm
Interior dimensions	w x d x h =* 53 x 45 x 128 cm	w x d x h =* 60 x 57 x 128 cm
Exterior dimensions with door open at 90°	w x d = 67 x 124 cm	w x d = 74 x 142 cm
Shelf size	w x d = 52,7 x 39 cm	w x d = 59 x 50 cm
Drawer inner dimensions	w x d x h = 50,4 x 32 x 5.6 cm	–
Max. load drawer/ shelf	16 kg/40 kg	–/40 kg
Weight	Net weight 90 kg, gross weight 104 kg	Net weight 113 kg, gross weight 127 kg

\*usable depth 5 cm less.

## Optional equipment

- **Glass door** with lock.
- **LED-Illumination** mounted on side wall.
- **Decorative door frame** for attaching the customer's decorative plates.
- **Door fitting.**
- **Equipped with castors.**
- **Condensate container** can be emptied manually if no condensation is required (e.g., in an operating theatre).
- **Plastic drawer on roller runners**
- **Additional length and cross dividers.**
- **Additional shelf** on rails or overlays.
- **Basket** on rails (only for LABO 468 PRO-ACTIVE) w x d x h = 60 x 45 x 10 cm.
- **Aluminium tray** with rails (only for LABO 468 PRO-ACTIVE) w x d x h = 60 x 40 x 2.5 cm.
- **Vario Support**, Storage Shelf on roller runners.
- **GSM-Module** for alarm text message transmission.
- **Temperature documentation:**
  - PC-KIT-NET (network capable version)
  - 4-pin PT 100 or PT 1000 temperature sensor, class 1/3 B, including built-in reference body (measurement bottle or refrigerating block).
  - Mechanical pen-recording thermometer with two waxed paper strips.